


Microcirculation: A Barometer of Tissue Health

Need: Quantitative assessment of tissue health



- Clinicians lack an on objective measure of circulation and tissue health at the point of care.
- **CURRENT STANDARD:** Subjective visual check, often by a non-expert.
- **NEED:** Rapid assessment of tissue health at the point-of-care



Solution: Spatial Frequency Domain Imaging (SFDI)



- Patented method to characterize tissue health
- Rapid (<1s) acquisition
- Non-contact and widefield
- Optical (non-radiative)
- Quantitative biomarkers



Solution Ecosystem



Clarifi® Imaging System

Cart-based, FDA-cleared microvascular imaging hardware outputs biomarkers to help identify circulatory issues below skin.

Magnifi Viewer

HIPAA-compliant cloud-based workflow to facilitate remote patient management by specialist care team.


Burn Digital Assessment

MTEC-22-08-BDA: Field-deployable device builds on established burn research foundation

- **Military need:** Handheld assessment tool to help healthcare teams triage and manage burn patients at far-forward point of care settings




1. Unique SFDI biomarkers correlate with histological burn depth at day 1



2. SFDI biomarkers enhance AI algorithms for predictive insights

Day 1 SFDI-based AI Output correlates with eventual outcome



- Unburned Skin
- Hyper-perfused/Superficial Burn
- Partial Thickness Burn
- Full Thickness Burn

Day 1 Multispectral AI Output did not correlate as closely with outcome



Year 1: Deliver a handheld prototype



- ~80% reduction in size¹
- ~70% reduction in weight¹
- ~5x reduction in cost of goods¹
- On-board processing
- On-board visualization

¹ compared to cart system imaging head


	Clarifi Cart System	m-Clarifi Handheld System
Imaging Head (weight)	5 lbs	1.5 lbs
Imaging Head (size)	20 cm x 19 cm x 10 cm	8.5 cm x 22.5 cm x 3.5 cm
Total System (weight)	105 lbs	1.5 lbs
Total System (size)	46 cm x 46 cm x 127 cm	8.5 cm x 22.5 cm x 3.5 cm
Processor	Surface tablet	On board compute platform
Camera	COTS	Custom
Light Source	COTS w/ motorized projector	Custom w/ no moving parts
Methods	Multi-image SFDI	Snapshot SFDI

Year 2+: Validation + clinical data collection

- Hardware Verification and Validation (2024)
- Software development (2024-2025)
- Clinical burn study at Institute of Surgical Research (2024-2025)
- FDA submission for broad clearance of biomarkers (2025)

Biomarkers: SFDI Platform Impact

Functional biomarkers have broad utility across the entire continuum of burn care



Outside DoD: Preventing limb complications

~700 Patients assessed



MICROVASCULAR ASSESSMENT

Patients w/potential circulatory issues identified remotely and referred to specialist

- Technician-acquired data in diabetes patients at clinical point-of-care
- Specialist remotely reviews cloud data to determine care path
- High risk patients are provided specialists care based on insights
- Proactive specialist care provided at the appropriate time & known to reduce complications by up to 50%

Limb care collaborators & publications

- ¹Lee, et al. *Journal of Diabetes and its Complications* (2020): 107624.
- ²Weinkauf, *Journal of Vascular Surgery* 69.2 (2019): 555-562.
- ³Murphy et al. *BMJ Open Diabetes Research and Care* (2020)
- ⁴Jett et al. *Journal of Diabetes Science and Technology* (2021)

